



Pre-shipment Inspection (PSI)

BASIC INSPECTION INFORMATION

Product:	400ml 3D dinosaur 3 in 1 shower gel – Red & Green & Orange	
Brand:	N/A	
Inspection location:	Qian xi industrial area, Lian xia town, Cheng hai district, Shan tou city,Guangdong province, China	
Inspection date:	Jul.08.2024	
Supplier:	Shan tou Kingkey industrial co.ltd	
Factory:	Shan tou Kingkey industrial co.ltd	
Inspector :	Colin tang	
Order quantity:	123552 PCS (13728 ctns)	
Qty ready for inspection (%):	46800PCS (37.8% finished, and 37.8% packed), (See attached qty breakdown)	
Reference sample available:	No	

INSPECTION RESULT JUDGMENT

OVERALL RESULT: Pending, see remark 1

CRITERIA	RESULTS	REMARKS			
1. Quantity	Pending	See remark 1			
2. Product and packing check	Passed	N/A			
3. On-site tests & special requirements	Passed	N/A			
4. Workmanship AQL:	Allowed quality level (AQL)		CRI	MAJ	MIN
		Found	0	11	11
		Max	0	14	21
		Checked	500	500	500

REMARKS

1.When leaving factory, the factory only finished and packed 46800 pcs (5000 ctns) for inspection.



INSPECTION STANDARDS	
Inspection type	Final Random Inspection (FRI)- <u>1st inspection</u>
Sampling standard	ANSI/ASQ Z1.4-2008 (MIL-STD-105E)
Sampling Plan (Single)	G-II (<u>normal inspection</u>)
AQL for Critical/ Major/ Minor	0 /1.5/4.0
Sampling size for Cri/Maj/Min	<u>500 /500 /500</u>
No. of cartons selected and inspected	No carton number printed on outer carton

1. Quantity	Result : Pending					
Approximate 80%of the total order quantity must be finished and packed	46800pcs(37.8 %) units finished and 46800_pcs (37.8%)packed into outer carton.					
	Item#	Oder required	Actual finished	Gift box packed	Outer carton packed	S/S
	Red dinosaur	41184	15600	15600	15600	166
	Green dinosaur	41184	15600	15600	15600	167
	Orange dinosaur	41184	15600	15600	15600	167
	Total	123552	46800	46800	46800	500

2.Product and packing check (photo)	Result : Passed
	
Finished products in warehouse	Finished products in warehouse



Outer carton view



Outer carton view- front & back



Outer carton view-sides



Outer carton packing with cardboard



Outer carton packing



PDQ packing view



PDQ packing view




PDQ packing view



PDQ packing view



PDQ packing view



Barcode sticker on PDQ



Scanning test



Marking on PDQ



Marking on PDQ



PDQ packing



Product view



Product view



Product view



Product view



3 in 1 sticker for orange dinosaur



3 in 1 sticker for red dinosaur



3 in 1 sticker for green dinosaur



Barcode sticker for orange dinosaur



Scanning test



Barcode sticker for red dinosaur



Scanning test



Barcode sticker for green dinosaur




Scanning test



Multiple layer stickers on bottom



Multiple layer stickers on bottom



Multiple layer stickers on bottom



Multiple layer stickers on bottom



Multiple layer stickers on bottom



Multiple layer stickers on bottom



Lip view



Lip view



Lip view



Lip view



Lip view



Lip view



Liquid valume check



Liquid valume check



Liquid valume check



Liquid valume check



Liquid valume check



Liquid valume check



Outer carton gross weight check




Outer carton length check



Outer carton width check




Outer carton height check



PDQ gross weight check



PDQ gross length check



PDQ gross width check



PDQ gross height check



Liquid color and scent check



Liquid color and scent check



Liquid color and scent check



Liquid bubble effect check



Liquid bubble effect check



Liquid bubble effect check



Liquid bubble effect check



Liquid bubble effect check



Liquid bubble effect check


Nil

N/A



Test Requirements	Comments from Inspector (Actual finding/photo)	Result
Universal test requirements		
1.Normal function check.(G-II: <u>500</u> pcs) (Check each function of the bottle lid, and all the function should be no problem)		Pass



2.Bar code scanning test (S-3: <u>32</u> pcs))	 	Pass
3. 3M tape adhesive test (S-3 <u>0</u> pcs) (Use 3M tape to do this test)	N/A	N/A
4. Internal structure/Major Components check (1-3 pcs per item)	N/A	N/A



5. Carton dimension and weight check (3 ctns per item) (The tolerance should be within +/- 5% of the carton size and weight which printed on carton)	Actual pallet size and weight measurement as below: <table border="1"> <thead> <tr> <th>Item#</th><th>Printed on master carton/Spec</th><th>Actual measuring</th></tr> </thead> <tbody> <tr> <td>Nil</td><td>36.1x26.9x18.0cm/ Nil</td><td>36.5x26.5x18.0cm/ 5.125kgs</td></tr> </tbody> </table>	Item#	Printed on master carton/Spec	Actual measuring	Nil	36.1x26.9x18.0cm/ Nil	36.5x26.5x18.0cm/ 5.125kgs	Pass
Item#	Printed on master carton/Spec	Actual measuring						
Nil	36.1x26.9x18.0cm/ Nil	36.5x26.5x18.0cm/ 5.125kgs						
6. Carton drop test (1 ctn per item) (ISTA-1A, 1 corner, 3 edge and 6 face of master carton, height : ____cm) , no master carton and product damaged.	The factory refused to do carton drop test during inspection due to the PDQ packing	N/A						
Special test requirements (Depend on different product)								



7. Transportation test for outer carton (1-2 ctns per item)

Use the factory's transportation test machine to do this test, and check if the product or package has any problem after test



Pass



8. Liquid leakage test (3 pcs per item)

Turn down the bottles for 2 hours , and then check if the Liquid is leaking or not



Pass



9. Bottle and liquid volume check (3 pcs per item)

Use the factory's volume test cup to check if the bottle volume is 400 ml, and then weight the liquid weight)



Pass



10. Liquid bubble effect check (3 pcs per item)

Mixing the shampoo with water, and then check the bubble effect



Pass



11. Smell test (3 pcs per item)
Compare the liquid sample, and check
scent of each item)



Pass



4. DEFECTIVES LIST

Defect description	CRI	MAJ	MIN	Defective Picture
Item# 400ml 3D dinosaur 3 in 1 shower gel – Red & Green & Orange				

Broken shrink film –热缩膜烂	0	9	0	
--------------------------	---	---	---	--







Scratching on the shrink film- 热缩膜刮伤	0	2	0	 
---	---	---	---	---



White mark(removable) -发白（可擦除）	0	0	5	  
-----------------------------------	---	---	---	---





Scratch mark –刮伤	0	0	4	 
------------------	---	---	---	--



Glue mark –胶水

0

0

1





Slanted 3 in 1 sticker – 3 合 1 贴纸歪	0	0	1	
---------------------------------------	---	---	---	--

Item# 400ml 3D dinosaur 3 in 1 shower gel – Red & Green & Orange				
Sample size	Maximum Allowed	Total Defects	Result	Notes
CRI : 500	0	0	Pass	N/A
MAJ :500	14	11	Pass	N/A
MIN :500	21	11	Pass	N/A

Remark: If no found the relevant AQL for the item reference , pls refer the AQL of most proximate sample size. for example, no AQL for sample size 150 pcs, pls refer the AQL of most proximate sample size of 125 pcs. If the AQL failed, the overall inspection result should be rejected, and pls make a remark like "Rejected due to AQL failed for item XXXX" on the fist page of report.

End of the Report