



Pre-shipment Inspection (PSI)

BASIC INSPECTION INFORMATION

Product:	Scent booster for laundry blue & pink	
Brand:	N/A	
Inspection location:	Guangdong province, China	
Inspection date:	Jun.07.2025	
Supplier:	GlamieaDailyGoods(Huizhou)Co.,Ltd	
Factory:	GlamieaDailyGoods(Huizhou)Co.,Ltd	
Inspector :	Colin tang & Lucky pan	
Order quantity:	79464 PCS (6622 ctns)	
Qty ready for inspection (%):	60864 PCS (76.5% finished, and 76.5 % packed), (See attached qty breakdown)	
Reference sample available:	Yes	

INSPECTION RESULT JUDGMENT

OVERALL RESULT: Pending, See remark 2,3

CRITERIA	RESULTS		REMARKS		
1. Quantity	Pending		See remark 2		
2. Product and packing check	Pass		N/A		
3. On-site tests & special requirements	Pending		See remark 3		
4. Workmanship AQL:	Allowed quality level (AQL)		CRI	MAJ	MIN
		Found	0	0	13
		Max	0	14	21
		Checked	500	500	500

REMARKS

1. During inspection, No customer's final approval sample for inspection reference.
2. During inspection, there are only 60864 pcs units (76.5%) were finished and packed for inspection, the factory promised to finish and pack all products before shipment.
3. During inspection, the outer carton transportation test can not be conducted due to their transportation test machine can not work .



INSPECTION STANDARDS							
Inspection type		Final Random Inspection (FRI)- <u>1st inspection</u>					
Sampling standard		ANSI/ASQ Z1.4-2008 (MIL-STD-105E)					
Sampling Plan (Single)		G-II (<u>normal inspection</u>)					
AQL for Critical/ Major/ Minor		0 /1.5/4.0					
Sampling size for Cri/Maj/Min		<u>500 /500 /500</u>					
No. of cartons selected and inspected		Pink 1481,1562,1571,1393,1540,1526,1510,3045,1308,3171,1457,3014,1076,2916,2772,1078,467,1574,1489,1535,1529 Blue 657,2076,386,428,254,647,245,128,770,1429,1523,352,2739,2684,546,402,3039,3022,2373,656,2073,2989					
1. Quantity		Result : Pending, see remark 2					
Approximate 80%of the total order quantity must be finished and packed		<u>60864</u> pcs(76.5%) units finished and <u>60864</u> pcs (76.5%)packed into outer carton.					
		Item#	Oder require d	Actual finished	Gift box packed	Outer carton packed	S/S
		Pink scent booster	39732	21840	21840	21840	250
		Blue scent booster	39732	39024	39024	39024	250
		Total	79464	60864	60864	60864	500
2.Product and packing check (photo)			Result : pass				
Item# Pink & Blue scent booster							
							
Finished products in wareho use			Outer carton view				



Outer carton view



Outer carton view



Outer carton packing



Outer carton packing



Gift box view



Gift box view



Gift box view



Gift box view




Gift box view



Gift box view



Gift box view



Gift box view



Gift box view





Gift box view




Gift box view



Barcode on gift box

扫描图片  条码类型 EAN_13 扫描时间 2025-06-07 15:39:18 扫描内容 复制内容 5903957310251	
Scanning test	Gift box view
	
Gift box view	Gift box view
	
Gift box view	Barcode on gift box
扫描图片  条码类型 EAN_13 扫描时间 2025-06-07 15:40:14 扫描内容 复制内容 5903957310244	
Scanning test	Gift box packing



Product view



Outer carton size check



Outer carton size check



Outer carton size check



Outer carton gross weight check



Product size check



Product size check



Product size check

3. ON-SITE TESTS & SPECIAL REQUIREMENTS

Result : Pending, see remark 3

Test Requirements	Comments from Inspector (Actual finding/photo)	Result
Universal test requirements		
1.Normal function check.(5_pcs_per item) (Put the scent booster inside the 25-40°C water, and then the scent booster must be dissolved within 20 Minutes, then use a white cotton fabric to test if the color is easy to dye on the fabric .)	 	Pass
2.Bar code scanning test (S-3: 32pcs))	 	Pass



3. 3M tape adhesive test (S-3 <u>0</u> pcs) (Use 3M tape to do this test)	N/A	N/A						
4. Internal structure/Major Components check (1-3 pcs per item)	N/A	N/A						
5. Carton dimension and weight check (3 ctns per item) (The tolerance should be within +/- 5% of the carton size and weight which printed on carton)	Actual pallet size and weight measurement as below: <table border="1"> <thead> <tr> <th>Item#</th><th>Printed on master carton/Spec</th><th>Actual measuring</th></tr> </thead> <tbody> <tr> <td>Pink & blue Scent booster</td><td>27.5x21x19cm/ 3.5 kgs</td><td>27.7x20.9x19cm/ 3.2kgs</td></tr> </tbody> </table>  	Item#	Printed on master carton/Spec	Actual measuring	Pink & blue Scent booster	27.5x21x19cm/ 3.5 kgs	27.7x20.9x19cm/ 3.2kgs	Pass
Item#	Printed on master carton/Spec	Actual measuring						
Pink & blue Scent booster	27.5x21x19cm/ 3.5 kgs	27.7x20.9x19cm/ 3.2kgs						



6. Carton drop test (1 ctn per item)
(ISTA-1A, 1 corner, 3 edge and 6 face of master carton, height : 76cm) , no master carton and product damaged.



Pass

Special test requirements (Depend on different product)

7. Transportation test for outer carton (1-2 ctns per item)
Use the factory's transportation test machine to do this test for 1 hour, and then check if the product or package has any problem after test

The factory's transportation test machine can not work during inspection

Pending

8. Product weight check (5 box per item)
The actual measured unit weight can not less than 180 g


Actual measured product weight as below (Unit:g):

Item#	1	2	3	4	5
Pink scent booster	182.8	183.8	187.3	184.5	185.2
	g	g	g	g	g
Blue scent booster	181.2	181.3	182.0	181.0	181.3
	g	g	g	g	g




Pass




9.Smell test (3 pcs per item) Compare the last approval sample, and check if the scent of each item is correct or not)		Pass
10. Dyeing test (1 pcs per item) Immerse a cotton fabric into the water with scent booster, and then wash the fabric in the clear water, then check if the color is stayed on the white fabric)		Pass



4. DEFECTIVES LIST

Defect description	CRI	MAJ	MIN	Defective Picture
Item# Pink scent booster				
There is a foreign body in the label-标签内有异物	0	0	1	



Round sticker is not centered-圆贴不居中	0	0	1	
Wrinkled round sticker-圆贴不平整	0	0	1	
Black stains-黑色污渍	0	0	1	
Item# Blue scent booster				



The round sticker is not centered-圆贴不居中	0	0	7	
--	----------	----------	----------	--



Batch number sticker is not centered-批号标签不居中	0	0	1	
edges of the round sticker is upturned-圆贴边缘翘起	0	0	1	

Item# Pink scent booster

Sample size	Maximum Allowed	Total Defects	Result	Notes
CRI : 250	0	0	Pass	N/A
MAJ :250	7	0	Pass	N/A
MIN : 250	14	4	Pass	N/A

Item# Blue scent booster

Sample size	Maximum Allowed	Total Defects	Result	Notes
CRI : 250	0	0	Pass	N/A
MAJ :250	7	0	Pass	N/A
MIN : 250	14	9	Pass	N/A

End of the Report

Remark: If no found the relevant AQL for the item reference , pls refer the AQL of most proximate sample size. for example, no AQL for sample size 150 pcs, pls refer the AQL of most proximate sample size of 125 pcs. If the AQL



failed, the overall inspection result should be rejected, and pls make a remark like "Rejected due to AQL failed for item XXXX" on the first page of report.